

FADS2 antibody

Product Information


Catalog No.:	FNab10288
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Fatty acid desaturase 2
Synonyms:	Acyl-CoA 6-desaturase Delta(6) fatty acid desaturase (D6D, Delta(6) desaturase, Delta-6 desaturase) Fatty acid desaturase 2 FADS2
Observed MW:	52 kDa
Uniprot ID :	O95864

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab10288(FADS2 antibody) at dilution of 1:1000