

## WDR74 antibody

### Product Information


Catalog No.:	FNab10301
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	WD repeat domain 74
Synonyms:	WD repeat-containing protein 74 NOP seven-associated protein 1 WDR74 NSA1
Observed MW:	40 kDa
Uniprot ID :	Q6RFH5

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200 - 1:1000
Image:	



K562 cells were subjected to SDS PAGE followed by western blot with FNab10301(WDR74 antibody) at dilution of 1:500