

CCL25 antibody

Product Information

Catalog No.:	FNab10310
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	chemokine(C-C motif) ligand 25
Synonyms:	C-C motif chemokine 25 Chemokine TECK Small-inducible cytokine A25 Thymus-expressed chemokine CCL25 SCYA25 TECK
Observed MW:	None
Uniprot ID :	O15444

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:200