

GYS2 antibody

Product Information


Catalog No.:	FNab10323
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	glycogen synthase 2
Synonyms:	Glycogen [starch] synthase, liver Glycogen synthase 2 GYS2
Observed MW:	81 kDa
Uniprot ID :	P54840

Application

Reactivity:	Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200 - 1:2000
Image:	



mouse liver were subjected to SDS PAGE followed by western blot with FNab10323(GYS2 Antibody) at dilution of 1:1000