

## SIK2 antibody

### Product Information


Catalog No.:	FNab10336
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	Serine/threonine-protein kinase SIK2
Synonyms:	Serine/threonine-protein kinase SIK2 Qin-induced kinase Salt-inducible kinase 2 (SIK-2) Serine/threonine-protein kinase SNF1-like kinase 2 SIK2 KIAA0781 QIK SNF1LK2
Observed MW:	104 kDa
Uniprot ID :	Q9H0K1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded rat heart tissue slide using FNab10336(SIK2 antibody) at dilution of 1:100



Hela cells were subjected to SDS PAGE followed by western blot with FNab10336(SIK2 antibody) at dilution of 1:1000