

TNFRSF6B antibody

Product Information

Catalog No.:	FNab10344
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Decoy receptor that can neutralize the cytotoxic ligands TNFS14/LIGHT, TNFSF15 and TNFSF6/FASL. Protects against apoptosis.

Immunogen information

Immunogen:	tumor necrosis factor receptor superfamily, member 6b, decoy
Synonyms:	Tumor necrosis factor receptor superfamily member 6B Decoy receptor 3 (DcR3) Decoy receptor for Fas ligand M68 TNFRSF6B DCR3 TR6
Observed MW:	33 kDa
Uniprot ID :	O95407

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HT29 cells were subjected to SDS PAGE followed by western blot with FNab10344(TNFRSF6B antibody) at dilution of 1:1000