

CAMK2D antibody

Product Information


Catalog No.:	FNab10387
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	calcium/calmodulin-dependent protein kinase II delta
Synonyms:	Calcium/calmodulin-dependent protein kinase type II subunit delta (CaM kinase II subunit delta, CaMK-II subunit delta) CAMK2D CAMKD
Observed MW:	30-40 kDa
Uniprot ID :	Q13557

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:500; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab10387(CAMK2D antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab10387(CAMK2D antibody) at dilution of 1:500