

ADIPOR1 antibody

Product Information


Catalog No.:	FNab10402
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3C3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	adiponectin receptor 1
Synonyms:	Adiponectin receptor protein 1 Progestin and adipoQ receptor family member 1 Progestin and adipoQ receptor family member 1 ADIPOR1 PAQR1 TESBP1A
Observed MW:	43 kDa
Uniprot ID :	Q96A54

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500 - 1:3000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab10402(ADIPOR1 antibody) at dilution of 1:1000