

## GREM1 antibody

### Product Information


Catalog No.:	FNab10411
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	Gremlin-1
Synonyms:	Gremlin-1 Cell proliferation-inducing gene 2 protein Cysteine knot superfamily 1, BMP antagonist 1 DAN domain family member 2 Down-regulated in Mos-transformed cells protein Increased in high glucose protein 2 (IHG-2) GREM1 CKTSF1B1 DAND2 DRM
Observed MW:	36 kDa
Uniprot ID :	O60565

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon using FNab10411(GREM1 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab10411(GREM1 antibody) at dilution of 1:1000