

KRT7 antibody

Product Information


Catalog No.:	FNab10415
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	keratin 7
Synonyms:	Keratin, type II cytoskeletal 7 Cytokeratin-7 (CK-7) Keratin-7 (K7) Sarcolectin Type-II keratin Kb7 KRT7 SCL
Observed MW:	54 kDa
Uniprot ID :	P08729

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab10415(KRT7 Antibody) at dilution of 1:50



Hela cells were subjected to SDS PAGE followed by western blot with FNab02217(KRT7 antibody) at dilution of 1:1000