

CAV1 antibody

Product Information


Catalog No.:	FNab10420
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	caveolin 1, caveolae protein, 22kDa
Synonyms:	Caveolin-1 CAV1 CAV
Observed MW:	50 kDa
Uniprot ID :	Q03135

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10420(Caveolin-1 Antibody) at dilution of 1:200



Various lysates were subjected to SDS PAGE followed by western blot with FNab10420(CAV1 antibody) at dilution of 1:1000