

GLB1 antibody

Product Information

Catalog No.:	FNab10432
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	galactosidase, beta 1
Synonyms:	Beta-galactosidase Acid beta-galactosidase (Lactase) Elastin receptor 1 GLB1 ELNR1
Observed MW:	75 kDa, 85 kDa
Uniprot ID :	P16278

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab10432(Beta galactosidase antibody) at dilution of 1:1000