

GATA4 antibody

Product Information


Catalog No.:	FNab10433
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	GATA binding protein 4
Synonyms:	Transcription factor GATA-4 GATA-binding factor 4 GATA4
Observed MW:	45-55 kDa
Uniprot ID :	P43694

Application

Reactivity:	Human,Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10433(GATA4 antibody) at dilution of 1:500