

ADAM28 antibody

Product Information

Catalog No.:	FNab10453
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in the adhesive and proteolytic events that occur during lymphocyte emigration or may function in ectodomain shedding of lymphocyte surface target proteins, such as FASL and CD40L. May be involved in sperm maturation.

Immunogen information

Immunogen:	ADAM metallopeptidase domain 28
Synonyms:	Disintegrin and metalloproteinase domain-containing protein 28 (ADAM 28) Epididymal metalloproteinase-like, disintegrin-like, and cysteine-rich protein II (eMDC II) Metalloproteinase-like, disintegrin-like, and cysteine-rich protein L (MDC-L) ADAM28 ADAM23 MDCL
Observed MW:	35-90 kDa
Uniprot ID :	Q9UKQ2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab10453(ADAM28 antibody) at dilution of 1:800