

FOXO6 antibody

Product Information

Catalog No.:	FNab10457
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcriptional activator.

Immunogen information

Immunogen:	forkhead box protein O6
Synonyms:	Forkhead box protein O6 FOXO6
Observed MW:	60 kDa
Uniprot ID :	A8MYZ6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10457(FOXO6 antibody) at dilution of 1:500