

NFATC3 antibody

Product Information

Catalog No.:	FNab10459
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as a regulator of transcriptional activation. Plays a role in the inducible expression of cytokine genes in T-cells, especially in the induction of the IL-2.

Immunogen information


Immunogen:	nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 3
Synonyms:	Nuclear factor of activated T-cells, cytoplasmic 3 (NF-ATc3, NFATc3) NFATx T-cell transcription factor NFAT4 (NF-AT4, NF-AT4c) NFATC3 NFAT4
Observed MW:	130-170 kDa
Uniprot ID :	Q12968

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab10459(NFATC3 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab10459(NFATC3 Antibody) at dilution of 1:1000