

## BIRC8 antibody

### Product Information

Catalog No.:	FNab10464
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Protects against apoptosis mediated by BAX.

### Immunogen information


Immunogen:	baculoviral IAP repeat-containing 8
Synonyms:	Baculoviral IAP repeat-containing protein 8 Inhibitor of apoptosis-like protein 2 (IAP-like protein 2, ILP-2) Testis-specific inhibitor of apoptosis BIRC8 ILP2
Observed MW:	30-40 kDa
Uniprot ID :	Q96P09

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab10464(BIRC8 antibody) at dilution of 1:100



various lysates were subjected to SDS PAGE followed by western blot with FNab10464(BIRC8 antibody) at dilution of 1:1000