

C3AR1 antibody

Product Information

Catalog No.:	FNab10502
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	C3a anaphylatoxin chemotactic receptor
Synonyms:	C3a anaphylatoxin chemotactic receptor (C3AR, C3a-R) C3AR1 AZ3B C3R1 HNFAG09
Observed MW:	55 kDa
Uniprot ID :	Q16581

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Rat spleen tissue were subjected to SDS PAGE followed by western blot with FNab10502(C3AR1 antibody) at dilution of 1:1000