

WFDC2 antibody

Product Information

Catalog No.:	FNab10508
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3F9
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	WAP four-disulfide core domain 2
Synonyms:	WAP four-disulfide core domain protein 2 Epididymal secretory protein E4 Major epididymis-specific protein E4 Putative protease inhibitor WAP5 WFDC2 HE4 WAP5
Observed MW:	None
Uniprot ID :	Q14508

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:500-1:2000
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab10508(HE4 Antibody) at dilution of 1:1000