

TMEM119 antibody

Product Information


Catalog No.:	FNab10517
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	transmembrane protein 119
Synonyms:	Transmembrane protein 119 Osteoblast induction factor (OBIF) TMEM119
Observed MW:	45 kDa
Uniprot ID :	Q4V9L6

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:200-1:1000
Image:	



Immunohistochemistry of paraffin-embedded human brain slide using FNab10517(TMEM119 antibody) at dilution of 1:500



Various lysates were subjected to SDS PAGE followed by western blot with FNab10517(TMEM119 antibody) at dilution of 1:500



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab10517(TMEM119 antibody) at dilution of 1:500