

MTNR1A antibody

Product Information


Catalog No.:	FNab10523
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Melatonin receptor type 1A
Synonyms:	Melatonin receptor type 1A (Mel-1A-R, Mella receptor) MTNR1A
Observed MW:	39 kDa
Uniprot ID :	P48039

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse brain slide using FNab10523(MTNR1A antibody) at dilution of 1:100



Rat kidney were subjected to SDS PAGE followed by western blot with FNab10523(FNab10523 antibody) at dilution of 1:1000