

ERCC4 antibody

Product Information


Catalog No.:	FNab10526
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	DNA repair endonuclease XPF
Synonyms:	DNA repair endonuclease XPF DNA excision repair protein ERCC-4 DNA repair protein complementing XP-F cells Xeroderma pigmentosum group F-complementing protein ERCC4 ERCC11 XPF
Observed MW:	110-120 kDa, 55-60 kDa
Uniprot ID :	Q92889

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



U-251 MG cells were subjected to SDS PAGE followed by western blot with FNab10526(XPF antibody) at dilution of 1:1000