

## SLC7A11 antibody

### Product Information


Catalog No.:	FNab10533
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	solute carrier family 7, (cationic amino acid transporter, y+ system) member 11
Synonyms:	Cystine/glutamate transporter Amino acid transport system xc- Calcium channel blocker resistance protein CCBR1 Solute carrier family 7 member 11 xCT SLC7A11
Observed MW:	35 kDa, 55 kDa
Uniprot ID :	Q9UPY5

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Rat liver were subjected to SDS PAGE followed by western blot with FNab10533(SLC7A11 antibody) at dilution of 1:1000