

IPO7 antibody

Product Information


Catalog No.:	FNab10626
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	importin 7
Synonyms:	Importin-7 (Imp7) Ran-binding protein 7 (RanBP7) IPO7 RANBP7
Observed MW:	130 kDa
Uniprot ID :	O95373

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:WB :	1:500-1:2000
Image:	



Various cell lysates were subjected to SDS PAGE followed by western blot with FNab10626(IPO7 antibody) at dilution of 1:1000