

TriMethyl-Histone H3-K27 antibody

Product Information


Catalog No.:	FNab10628
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	A synthetic trimethylated peptide around K27 of human histone H3
Synonyms:	Histone H3.1t (H3/t, H3t) H3/g Histone H3.4 H3-4 H3FT HIST3H3
Observed MW:	16 kDa
Uniprot ID :	Q16695

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:	WB : 1:500-1:1000
Image:	



Various nuclear extract were subjected to SDS PAGE followed by western blot with FNab10628 at dilution of 1:800