

## GDF15 antibody

### Product Information


Catalog No.:	FNab10629
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	growth differentiation factor 15
Synonyms:	Growth/differentiation factor 15 (GDF-15) Macrophage inhibitory cytokine 1 (MIC-1) NSAID-activated gene 1 protein (NAG-1) NSAID-regulated gene 1 protein (NRG-1) Placental TGF-beta Placental bone morphogenetic protein Prostate differentiation factor GDF15 MIC1 PDF PLAB PTGFB
Observed MW:	34 kDa
Uniprot ID :	Q99988

### Application

Reactivity:	Human
Tested Application:	WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta using FNab10629(GDF15 antibody) at dilution of 1:100



HT1080 cell lysates were subjected to SDS PAGE followed by western blot with FNab10629(GDF15 antibody) at dilution of 1:1000