

HBc antibody

Product Information

Catalog No.:	FNab10631
Size:	1mg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7H1
IsoType:	IgG1
Storage:	PBS with 0.1% Proclin300 pH 7.3, store at 2-8°C up to 21 days, and at -20°C for long-term storage. (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	hepatitis B virus core antigen
Synonyms:	None
Observed MW:	None
Uniprot ID :	None

Application

Reactivity:	None
Tested Application:	Recommended for ELISA, RIA and CLIA. Each laboratory should determine an optimum working titer for use in its particular platform application.