

## **RAD50 antibody**

### **Product Information**


Catalog No.:	FNab10634
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### **Immunogen information**


Immunogen:	RAD50 homolog
Synonyms:	DNA repair protein RAD50  hRad50  RAD50  RAD50 2  RAD50 homolog
Observed MW:	160 kDa
Uniprot ID :	Q92878

### **Application**

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab10634(RAD50 Antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab10634(RAD50 Antibody) at dilution of 1:1000