

NEU2 antibody

Product Information


Catalog No.:	FNab10637
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	N-acetyl-alpha-neuraminidase 2
Synonyms:	Sialidase-2 Cytosolic sialidase NEU2 N-acetyl-alpha-neuraminidase 2
Observed MW:	42 kDa
Uniprot ID :	Q9Y3R4

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Raji cells were subjected to SDS PAGE followed by western blot with FNab10637(NEU2 Antibody) at dilution of 1:1000