

BRWD1 antibody

Product Information


Catalog No.:	FNab10640
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Bromodomain and WD repeat-containing protein 1
Synonyms:	BRWD1 C21orf107 FLJ43918 N143 WD repeat containing protein 9 WDR9
Observed MW:	263 kDa
Uniprot ID :	Q9NSI6

Application

Reactivity:	Human, Mouse
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10640(BRWD1 Antibody) at dilution of 1:1000