
RAD54B antibody

Product Information


Catalog No.:	FNab10644
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	DNA repair and recombination protein RAD54B
Synonyms:	DNA repair and recombination protein RAD54B, RAD54 homolog B, RAD54B
Observed MW:	103 kDa
Uniprot ID :	Q9Y620

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10644(RAD54B Antibody) at dilution of 1:1000