

---

## HRPT2 antibody

### Product Information


Catalog No.:	FNab10645
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7C9
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	cell division cycle 73, Paf1/RNA polymerase II complex component, homolog
Synonyms:	C1orf28, CDC73, HPT JT, HRPT2, HRPT2, CDC73, Hyperparathyroidism 2 protein, Parafibromin
Observed MW:	61 kDa
Uniprot ID :	Q6P1J9

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab10645(HRPT2 antibody) at dilution of 1:50



Various lysates were subjected to SDS PAGE followed by western blot with FNab10645(HRPT2 antibody) at dilution of 1:1000