

---

## SLC12A9 antibody

### Product Information


Catalog No.:	FNab10651
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Solute carrier family 12 member 9
Synonyms:	CCC interacting protein 1, CCC6, CIP1, hCCC6, hCIP1, SLC12A9, WO3.3
Observed MW:	100 kDa
Uniprot ID :	Q9BXP2

### Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:WB:	1:500-1:1000
Image:	



HEK293 cells were subjected to SDS PAGE followed by western blot with FNab10651(SLC12A9 Antibody) at dilution of 1:500