
OPA1 antibody

Product Information


Catalog No.:	FNab10658
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Optic atrophy protein 1
Synonyms:	KIAA0567, largeG, MGM1, NTG, OPA1, Optic atrophy protein 1
Observed MW:	80-100 kDa
Uniprot ID :	O60313

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10658(OPA1 antibody) at dilution of 1:1000