

---

## JAK3 antibody

### Product Information


Catalog No.:	FNab10659
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Tyrosine-protein kinase JAK3
Synonyms:	JAK 3, JAK3, JAK3_HUMAN, JAKL, Janus kinase 3, L JAK, Leukocyte janus kinase, LJAK, Tyrosine protein kinase JAK3
Observed MW:	125 kDa
Uniprot ID :	P52333

### Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HEK293 cells were subjected to SDS PAGE followed by western blot with FNab10659(JAK3 antibody) at dilution of 1:1000