

---

## TREM2 antibody

### Product Information


Catalog No.:	FNab10660
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Triggering receptor expressed on myeloid cells 2
Synonyms:	TREM 2, TREM2, Trem2a, Trem2b, Trem2c
Observed MW:	40 kDa
Uniprot ID :	Q9NZC2

### Application

Reactivity:	Human, Mouse
Tested Application:	WB
Recommended dilution: WB:	1:500-1:2000
Image:	



THP-1 cells were subjected to SDS PAGE followed by western blot with FNab10660(TREM2 antibody) at dilution of 1:1000