
PPARG antibody

Product Information


Catalog No.:	FNab10662
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2C10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	peroxisome proliferator-activated receptor gamma
Synonyms:	CIMT1, NR1C3, PPAR gamma, PPARG, PPARG1, PPARG2, PPARgamma, PPAR γ
Observed MW:	50-60 kDa
Uniprot ID :	P37231

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab10662(PPARG Antibody) at dilution of 1:100



various lysates were subjected to SDS PAGE followed by western blot with FNab10662(PPARG antibody) at dilution of 1:1000