

NLRP3 antibody

Product Information


Catalog No.:	FNab10665
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3H9
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	NLR family, pyrin domain containing 3
Synonyms:	AGTAVPRL, AII, AII/AVP, AVP, C1orf7, Caterpillar protein 1.1, CIAS1, CLR1.1, Cryopyrin, FCAS, FCU, MWS, NALP3, NLRP3, PYPAF1
Observed MW:	118 kDa
Uniprot ID :	Q96P20

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10665(NLRP3 antibody) at dilution of 1:1000