

---

## GRIN1 antibody

### Product Information


Catalog No.:	FNab10666
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1J4
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Glutamate receptor ionotropic, NMDA 1
Synonyms:	GRIN1, GluN1, MRD8, NMD-R1, NMDA1, NMDAR1, NR1, NMDA 1
Observed MW:	120 kDa
Uniprot ID :	Q05586

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:5000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10666(GRIN1 antibody) at dilution of 1:2000