

---

## XRCC1 antibody

### Product Information


Catalog No.:	FNab10668
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	X-ray repair complementing defective repair in Chinese hamster cells 1
Synonyms:	DNA repair protein XRCC1, RCC, XRCC1
Observed MW:	80-90 kDa
Uniprot ID :	Q13315

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded rat testis tissue slide using FNab010668( XRCC1 antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab10668(XRCC1 antibody) at dilution of 1:1000