
Phospho-PERK(Thr982) antibody

Product Information


Catalog No.:	FNab10673
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	eukaryotic translation initiation factor 2-alpha kinase 3
Synonyms:	EIF2AK3, HRI, HsPEK, Pancreatic eIF2 alpha kinase, PEK, PERK, WRS
Observed MW:	180 kDa
Uniprot ID :	Q9NZJ5

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:1000-1:5000
Image:	



HEK293 cells were subjected to SDS PAGE followed by western blot with FNab10673 at dilution of 1:2000