
HNRNPM antibody

Product Information


Catalog No.:	FNab10674
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Heterogeneous nuclear ribonucleoprotein M
Synonyms:	HNRNPM, HNRNPM4, HNRPM, HNRPM4, HTGR1, NAGR1
Observed MW:	73-78 kDa
Uniprot ID :	P52272

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab10674(HNRNPM antibody) at dilution of 1:1000