
CLH1 antibody

Product Information


Catalog No.:	FNab10675
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Clathrin heavy chain 1
Synonyms:	CHC, CHC17, Clathrin heavy chain 1, clathrin, heavy chain (Hc), CLH17, CLH17, CLTC, CLTCL2, Hc, KIAA0034
Observed MW:	180 kDa
Uniprot ID :	P52272

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10675 (CLH1 antibody) at dilution of 1:1000