

---

## SLC30A1 antibody

### Product Information


Catalog No.:	FNab10680
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Proton-coupled zinc antiporter SLC30A1
Synonyms:	SLC30A1, Zinc transporter 1, ZnT 1, ZNT1, ZRC1
Observed MW:	60-70 kDa
Uniprot ID :	Q9Y6M5

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



THP-1 cells were subjected to SDS PAGE followed by western blot with FNab10680(SLC30A1 antibody) at dilution of 1:1000