

---

## IL1RAP antibody

### Product Information


Catalog No.:	FNab10686
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Interleukin-1 receptor accessory protein
Synonyms:	C3orf13  FLJ37788  IL 1R 3  IL 1R3  IL 1RAcP  IL1R3  IL1RAP  Interleukin 1 receptor 3
Observed MW:	70 kDa
Uniprot ID :	Q9NPH3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



various lysates were subjected to SDS PAGE followed by western blot with FNab10686(IL1RAP antibody) at dilution of 1:1000