
Acetyl-Histone H3(K27) antibody

Product Information


Catalog No.:	FNab10687
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	A synthetic peptide
Synonyms:	Histone H3.1 Histone H3/a Histone H3/b Histone H3/c Histone H3/d Histone H3/f Histone H3/h Histone H3/i Histone H3/j Histone H3/k Histone H3/l H3C1 H3FA HIST1H3A
Observed MW:	15 kDa
Uniprot ID :	P68431

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000-1:5000
Image:	



various lysates were subjected to SDS PAGE followed by western blot with FNab10687 at dilution of 1:2000