
PLA2G4A antibody

Product Information


Catalog No.:	FNab10689
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3J9
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Cytosolic phospholipase A2
Synonyms:	Cytosolic phospholipase A2 cPLA2 Phospholipase A2 group IVA Phospholipase A2 Phosphatidylcholine 2-acylhydrolase PLA2G4A
Observed MW:	85-100 kDa
Uniprot ID :	P47712

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000 - 1:5000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10689(PLA2G4A antibody) at dilution of 1:2000