

---

## GCH1 antibody

### Product Information


Catalog No.:	FNab10692
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	GTP cyclohydrolase 1
Synonyms:	DYT14  DYT5  GCH  GCH1  GTP CH 1  GTP CH I  GTP cyclohydrolase 1  GTP cyclohydrolase I  GTPCH1
Observed MW:	30 kDa
Uniprot ID :	P30793

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10692(GCH1 antibody) at dilution of 1:1000