

## CCNK antibody

### Product Information


Catalog No.:	FNab10694
Size:	100µg
Form:	liquid
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	CCNK CPR4 cyclin K
Uniprot ID :	O75909

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000



Hela lysates were subjected to SDS PAGE followed by western blot with FNab10694(CCNK antibody) at dilution of 1:1000