
D-Dimer antibody

Product Information

Catalog No.:	FNab10695
Size:	100µg/1mg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	5D8
IsoType:	IgG1
Storage:	10mMPBS pH 7.4, 0.1% proclin300, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Application

Reactivity:	Human
Tested Application:	CLIA, GICA(capture antibody)